

NASH-FibroTest

Estimating NASH, Steatosis and Fibrosis in NAFLD (non-alcoholic liver disease) patients

NASH-FibroTest is a liver panel including 3 non-invasive tests: NashTest 2, SteatoTest 2 and FibroTest, along with 2 additional tests: ActiTest and AshTest.

Rational for NASH-FibroTest (NashTest 2 + SteatoTest 2 + FibroTest)

NASH-FibroTest is a liver panel used in the diagnosis and the follow-up of liver NASH, steatosis and fibrosis/cirrhosis with one standard blood sample, done at a local laboratory:

- NashTest 2: estimates the liver inflammation as a quantitative assessment of steatohepatitis and the prediction of liver outcome²
- SteatoTest 2: estimates the liver steatosis³
- FibroTest: estimates the liver fibrosis⁵, predicts the long-term liver-related mortality⁶

Clinicians need the 3 diagnoses to take a decision for treatment and/or followup.⁷⁻⁸

NashTest 2

NashTest 2 is constructed using updated histological consensus on NASH definition.9 It is a quantitative test, assessing the severity of the liver inflammation (NASH), without the need of BMI.

The absence of glucose assay makes NashTest 2 more robust in patients with type 2 diabetes.

SteatoTest 2

Hepatic steatosis, which is assessed by the SteatoTest 2, is a build-up of fat in the liver, which frequently causes elevated levels of gamma-GT and transaminases.¹⁰ SteatoTest is recommended by EASL, EASD and EASO for the management of nonalcoholic liver disease.11 SteatoTest 2 is equivalent to the SteatoTest (noninferiority) but is more convenient (no BMI).3

FibroTest:

Reference test

FibroTest is recommended by the EASL, EASD and EASO for the management of nonalcoholic liver disease.11 FibroTest is the reference test for liver fibrosis and cirrhosis (recommended by AASLD10, EASL-ALEH¹², CASL¹³, APASL¹⁴ and WHO¹⁵ in chronic liver diseases).

Why monitor the liver of NAFLD patients?

NASH-FibroTest evaluates NASH¹, steatosis³, fibrosis (incl. cirrhosis)⁵ and therefore the risk of liver cancer in populations with either: metabolic-related factors (i.e. overweight, hypertension, dyslipidemia or type 2 diabetes)^{5,16-17} with coronary artery disease⁴ or with ultrasound steatosis⁵.

NASH-FibroTest comes with 2 additional tests


ActiTest and AshTest will also be provided along with NASH-FibroTest as additional tests since they may be very useful when NAFLD is associated with chronic viral hepatitis (ActiTest) or alcohol consumption (AshTest).

	FibroTest	Transient Elastography	APRI	FIB-4
Applicability				
All stages performance			•	
Inflammatory activity false positive	•	•	•	•
Observer variability		•	•	•
Cost		•		
Prognosis				

FibroTest benchmark in direct comparisons, in intention to diagnose 18

NASH-FibroTest is the only package on the market combining the EASL recommended non-invasive bloodbased tests for both steatosis and fibrosis.


A complete NASH panel based on a blood sample (no BMI)





NashTest 2 assesses NASH

- ▶ N0 : no NASH
- NI: mild NASH
- ▶ N2 : moderate NASH
- N3 : severe NASH




0.00



SteatoTest 2 assesses liver steatosis

- ▶ S0 : no steatosis (<5%)
- ▶ SI: mild steatosis (but clinically significant) (5-33%)
- ▶ S2S3 : moderate to severe steatosis (clinically significant) (34-100%)





FibroTest estimates liver fibrosis:

- ▶ F0: no fibrosis
- ▶ F1: minimal fibrosis
- ▶ F2: moderate fibrosis
- F3: advanced fibrosis
- ▶ F4: severe fibrosis (cirrhosis)

Additional tests:





to alcohol:

ActiTest estimates necroinflammatory activity:

- A0: no activity A1: minimal activity
- A2: moderate activity A3: severe activity
- H0 : no alcoholic steatohepatitis HI: moderate alcoholic steatohepatitis
- ▶ H2 : significant alcoholic steatohepatitis

AshTest assesses inflammation due

H3: severe alcoholic steatohepatitis

Assays (done at a local lab.): Alpha-2

macroglobulin, Haptoglobin, Apolipoprotein AI, Total Bilirubin, GGT, ALT, AST, Cholesterol, Triglycerides, Fasting Glucose, age, sex- according to BioPredictive precautions of use (biopredictive.com)

- Poynard T. et al. Eur | Gastroenterol Hepatol. 2018 (May)
- Poynard T. et al. Eur.] Gastroenterol Hepatol. 2018 (May) Munteanu M. et al. Aliment Pharmacol Ther. 2018 Poynard T. et al. Eur.] Gastroenterol Hepatol. 2018 (Dec.) Perazzo H. et al. Aliment Pharmacol Ther. 2016 Munteanu M. et al. Aliment Pharmacol Ther. 2016 Munteanu M. et al. Aliment Pharmacol Ther. 2018 Bedossa Pet al. Hepatology 2012 Bedossa Pet al. Hepatology 2014 Poynard T. et al. Eur.] Gastroenterol Hepatol. 2018 (April) Chalasani N. et al. Hepatology 2018 EASL-EASD-EASO Clinical Guidelines J. Hepatol. 2016 EASL-ALEH Panel. J. Hepatol. 2015 CASL Consensus guidelines Can.] Gastroenterol. 2012 Shiha G. et al. Hepatol Int. 2009 Guidelines WHO 2014 and 2015 Poynard T. et al. Plos One 2012 Jacquerninet S. et al. Clin Gastroenterol Hepatol. 2008 Houot M. et al. Aliment Pharmacol Ther. 2016

